


Phenion® FT Skin Model
RNA Isolation I
 (Mechanical homogenization)



Preface

Reliable results in RNA-preparation are routinely achievable when using RNeasy® Mini Kits from Qiagen (Hilden, Germany). For safety instructions, please refer the manufactures instructions. Please note that other suppliers might also lead to excellent results.

For maximum RNA-yields, we recommend to homogenize the tissue using a mixer mill and to follow this protocol – which is slightly adapted from the Qiagen RNeasy Mini Kit protocol. If a mixer mill is not available, we recommend using the alternative protocol “Phenion® FT Skin Model RNA-Isolation II”.

Typically, 25-40 µg RNA, with an OD_{260/280} ratio of ≥ 2.0 , is obtained from one half of the Phenion® FT Skin Model. This yield suits most RNA analyses methods, enabling to subject the other half of the tissue to further analytical methods, e.g. histological assessments or protein analyses.

Materials/Disposables/Reagents

Material	Company	Order-No.
RNeasy Mini Kit	Qiagen, (Hilden, Germany)	74104
Disposable Scalpel	e.g. Braun (Tuttlingen, Germany)	5518059
1.5 ml / 2.0 ml Reaction tubes (PCR-grade, RNase free)	Eppendorf (Hamburg, Germany)	0030 123.328 0030 123.344
Mixer Mill e.g. Tissue Lyzer	Qiagen, (Hilden, Germany)	85300
Centrifuge e.g. Biofuge fresco	Heraeus Kendro, (Osterode, Germany)	750 05521

Procedure

1. Dissect Phenion® FT Skin Model in two halves (store first half appropriately for further processing of your choice) and cut second half into eight equally sized pieces.
2. Transfer tissue pieces into 2 ml reaction tube and add 350 µl RLT buffer supplemented with 10 µl β-Mercaptoethanol/ml buffer. Insert 5 mm stainless steel beads and close vial.
3. Disrupt tissue with a Mixer Mill for 5 min with agitation speed of 30 Hz
4. Pipet homogenate into a new tube. Add 500 µl deionized water and 10 µl proteinase K. Mix gently and incubate 10 min at 55 °C.
5. Pellet cell debris by centrifugation at 8.000 g for 30 sec.
6. Pipet supernatant into a new collection tube and add 0.5 volumes 98% ethanol. Mix gently by pipetting up and down slowly (do not vortex!).
7. Transfer 700 µl onto RNeasy Mini Kit Spin Column. Centrifuge for 15 sec at 8.000 x g and discard flow-through.
8. Transfer remaining lysate onto corresponding spin column. Centrifuge for 15 sec at 8.000 x g and discard flow-through.
9. Wash spin column with 350 µl RW1 buffer. Centrifuge for 15 sec at 8.000 g and discard flow-through.
10. Add 80 µl DNase solution (10 µl DNase stock solution + 70 µl RDD buffer) to the spin column membrane and incubate 15 min at RT.
11. Add 350 µl RW1 buffer to the spin column. Centrifuge for 30 sec at 8.000 x g and discard flow-through.
12. Add 500 µl RPE buffer. Centrifuge for 30 sec at 8.000 x g and discard flow-through.
13. Place spin column on a new 2 ml collection tube and repeat washing step with 500 µl RPE buffer.
14. Dry membrane of spin column by centrifugation for 4 min at max speed.
15. Place spin column carefully in a new 1.5 ml collection tube and add 50 µl RNase-free water on the surface of the membrane. Centrifuge for 3 min at 8.000 x g to elute the RNA and store appropriately.