


E-Mail: info@phenion.com

Web: www.phenion.com

Phenion® Cell Culture Media Contents



Component	Fibroblast Growth Medium	Keratinocyte Growth Medium	Air-liquid Interface Growth Medium
DMEM with Glutamax	✓	✓	✓
Ham/F-12 Medium		✓	✓
L-Ascorbic acid 2-phosphate	✓	✓	✓
Penicillin/Streptomycin	✓	✓	✓
Hydrocortisone		✓	✓
Insulin		✓	✓
Bovine serum albumin (BSA)			✓
FetalClone II serum		✓	
Fetal calf serum	✓		
Cholera toxin		✓	
Epidermal growth factor		✓	
Triiodo-L-thyronine		✓	
Adenine		✓	

Ref.: K.R. Mewes et al., *Skin Pharmacol. Physiol.*, 2007, 20:85-95 – DOI: 10.1159/000097655